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**Nozzle**

100- $\mu$ m nozzle is removable and can be sonicated.

A registered key-fit position at the bottom of the cuvette provides fixed stream alignment.

**Fluidic tanks**

Autoclavable 10-L stainless steel sheath container

10-L polypropylene waste container

**Performance****System operation**

Time from system powerup to running sample, including powerup, fluidic startup, stream optimization and optional performance checks (instrument and drop delay QC): <20 minutes; excluding optional performance checks: <10 minutes

Sorting nozzle tip can be removed during operation, replaced and stream optimized: <3 minutes

**Fluorescence sensitivity**

Fluorescence sensitivity was measured using SPHERO™ Rainbow Calibration Particles according to the manufacturer's specifications:

FITC: <110 molecules of equivalent soluble fluorochrome (MESF-FITC)

PE: <30 molecules of equivalent soluble fluorochrome (MESF-PE)

**Fluorescence detection efficiency**

Qr is the relative fluorescence detection efficiency, used for describing the detection efficiency of a detector for a specific fluorochrome. The units for Qr are photoelectrons generated per assigned BD unit (ABD unit). One ABD unit, for a given fluorochrome, is defined as the fluorescence of one antibody (fluorochrome to protein 1:1) bound to a cell. The higher the Qr value, the better the relative fluorescence detection efficiency per molecule reference fluorochrome measured in that detector.

FITC Qr (x1,000): 40 photoelectrons/ABD

PE Qr (x1,000): 325 photoelectrons/ABD

CD4 FITC: 1,600 photoelectrons\*

CD4 PE: 13,000 photoelectrons\*

\*CD4 antibodies bound = 40,000

Qr values from the same lot of beads were taken from one BD FACSMelody Instrument. Qr values can vary between instruments and instrument configurations.

**Fluorescence resolution**

Full-peak coefficient of variation (FPCV): <3.0%, G0/G1 peak for propidium iodide (PI)-stained chicken erythrocyte nuclei (CEN)

**Fluorescence linearity**

Doublet/singlet ratio. PI-stained CEN: 1.95–2.05

**Forward and side scatter sensitivity**

Sensitivity enables separation of 0.5- $\mu$ m beads from noise.

**Forward and side scatter resolution**

Scatter performance is optimized for resolving lymphocytes, monocytes and granulocytes.

**Sort performance****Droplet sorting**

34,000 drops per second

Automated setup, optimization and monitoring of droplet breakoff and sort streams

Automated drop-delay determination with BD FACST™ Accudrop beads

Automated clog detection and sort tube protection system using Sweet Spot technology

**Purity and yield**

At 23 psi and 34 kHz with an average of 10,000 events per second, a 4-way sort of 5% target populations achieved a purity of 98% and a yield of >80% of Poisson's expected yield for all four populations. No limitation is imposed on the event rate. However, yield decreases with higher event rates based on Poisson distribution principles. Similar results are observed for 2-way sorts.

#### Sort collection

4-way sorting: 1.5-, 2.0-, and 5.0-mL tubes

2-way sorting: 1.5-, 2.0- and 5.0-mL tubes

1-way sorting: 6-, 24-, 48-, 96- and 384-well plates, 96-well PCR tray, microscope slide (optional)

Temperature control: water recirculation unit to provide heating or cooling for collection into tube holders, multiwell plates and slides (optional)

#### Index Sorting

Correlation of flow cytometry parameters of sorted events with well location of a multiwell plate, tray or slide (standard feature enabled with optional plate sorting).

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